

MANUAL HANDLING

HoverMatt® Single-Patient Use



HoverMatt® Single-Patient Use

The number one choice for lateral patient transfers now offers the added convenience of disposability in the Single-Patient Use HoverMatt® transfer mattress. The inflatable transfer mattress works just like the reusable HoverMatt® air transfer system, using a cushion of air that enables caregivers to safely transfer patients without lifting or straining. By eliminating injuries related to lateral transfers and repositioning, the Single-Patient Use HoverMatt® can help improve staff retention while meeting legislative guidelines for safe patient handling.



Features:

- Improves patient and caregiver safety and comfort during transfers
- Reduces skin shear and bruising
- Eliminates reprocessing (no laundering)
- Ideal for use in any patient care department
- Radiolucent/MRI compatible
- Available in multiple styles and widths to accommodate all patient sizes
- Significantly reduces the possibility of cross contamination
- Reduces costs and logistics associated with reprocessing
- Variable speed air supply provides rapid inflation
- Weight Limit 544kg

SHA010914-2.4

MANUAL HANDLING

HoverMatt® Single-Patient Use



Product Specifications

Material Specs

Sewn Construction

Top: Non-Woven Polypropylene Fibre

Bottom: Nylon Fabric

Latex free

Water Soluble Perimetre thread

(120°) - Thread in sewn handles

in NOT Water Soluble

Available in Three Sizes

86cm W x 198cm L - MHHT-HMSPU-34

99cm W x 198cm L - MHHT-HMSPU-39

127cm W x 198cm L - MHHT-HMSPU-50

86cm W x 115cm L - MHHT-HMSPU-34MI

99cm W x 115cm L - MHHT-HMSPU-39MI

MHHT-HMSPU-34, MHHT-HMSPU-39, MHHT-HMSPU-34MI, MHHT-HMSPU-39MI


10 per box

MHHT-HMSPU-50

5 per box

Compatible with HoverTech Air Supplies

MHHT-AIR2300 (Variable Speed - 230 V)



Product Code	Product Name
MHHT-HMSPU-34EA	HoverMatt Lateral Transfer Single Patient Use 34" - Each
MHHT-HMSPU-34	HoverMatt Lateral Transfer Single Patient Use 34" - Box of 10
MHHT-HMSPU-39EA	HoverMatt Lateral Transfer Single Patient Use 39" - Each
MHHT-HMSPU-39	HoverMatt Lateral Transfer Single Patient Use 39" - Box of 10
MHHT-HMSPU-50EA	HoverMatt Lateral Transfer Single Patient Use 50" - Each
MHHT-HMSPU-50	HoverMatt Lateral Transfer Single Patient Use 50" - Box of 5



SHA010914-2.4