

MANUAL HANDLING

HoverMatt® Mini (Half-Matt)



HoverMatt® Mini (Half-Matt)

The HoverMatt® Mini (Half-Matt) delivers the same effortless lateral transfers and positioning as the standard HoverMatt® Air Transfer System, but in a shorter length. Specifically designed for use in the O.R. with specialty tables as well as in Labour and Delivery, the Mini is available in both reusable and single-patient use models. The reusable Mini with a double-coated finish is often preferred when excessive fluids make stain resistance a priority.



Variable Speed
HoverMatt® Air Supply

Double Coated
HoverMatt® Mini (Half-Matt)

Single Patient Use
HoverMatt® Mini (Half-Matt)

Features:

- Air-assisted technology improves patient and caregiver safety and comfort during transfers
- Inflated mattress moves patient in stable position, reducing skin shear and bruising
- Half-length design is ideal for use with O.R. specialty tables
- Radiolucent and MRI compatible for artifact-free imaging
- Available in multiple styles and widths with a 544kg weight limit to accommodate most patients
- Variable speed air supply allows quieter and slower inflation options

SHA010914-2.2

MANUAL HANDLING

HoverMatt® Mini (Half-Matt)



Product Specifications

Single Patient Use Mini (Half Matt)

Sewn Construction
 Top: Non-woven polypropylene fiber
 Bottom: Nylon fabric
 Water soluble perimetre thread (120)
 Thread in sewn handles is NOT water soluble
 Latex free

Available Size:

86cm W x 114cm L - MHHT-HMSPU-34MI Box of 10
 99cm W x 114cm L - MHHT-HMSPU-39MI Box of 10

Heat Sealed Mini (Half Matt)

Heat-Sealed Construction: RF welding in all seams
 Material: Nylon twill
 Latex free

Available Size:

86cm W x 114cm L - MHHT-HMHS-34MI Each

Double Coated Mini (Half Matt)


Heat-Sealed Construction: RF welding in all seams
 Material: Nylon twill
 Latex free

Available Size:

86cm W x 114cm L - MHHT-HMDC-34MI Each

Compatible with HoverTech Air Supplies

MHHT-AIR2300 (Variable Speed - 230 V)



If a Single Patient Use Mini (Half Matt) transfer mattress is laundered, it will no longer inflate when connected to the air supply. The 'DO NOT WASH' label at the footend of the product will dissolve and state, 'DISCARD IMMEDIATELY'

Product Code	Product Name
MHHT-HMSPU-34MI	HoverMatt Lateral Transfer Single Patient Use 34" MINI - Box of 10
MHHT-HMSPU-39MI	HoverMatt Lateral Transfer Single Patient Use 39" MINI - Box of 10
MHHT-HMHS-34MI	HoverMatt Lateral Transfer Heat Sealed 34" MINI - Each
MHHT-HMDC-34MI	HoverMatt Lateral Transfer Double Coated 34" MINI - Each



SHA010914-2.2