

Proxi-Mate



www.statinahealthcare.com.au

Proxi-Mate Pro Wireless Set

The most recent addition to the range, **Proxi-Mate Pro Wireless Set** is a falls monitoring system specifically designed to meet the needs of a busy Acute or Sub Acute Rehabilitation Unit.

All-round use for monitoring residents in aged care facilities, neurologically compromised patients in acute care and people affected by dementia and other conditions, in their own homes. Can operate wireless.

The system runs on its own rechargeable batteries so it is not tied to any location. It will 'talk' to up to 8 different pagers and the 16 character message can be set to read whatever is needed.

The system records the last 4000 events, including time and date events occur.*

* 100 events viewable on screen, full log can be downloaded to a USB drive via the **Proxi-Log** device.

Proxi-Mate Pro Nurse Call Set

The tough, die cast aluminium **Proxi-Mate Pro Nurse Call Set** is designed to plug into most existing nurse call systems commonly in use today. The units are ideal for the Nursing Home environment where simplicity can make the difference.

A simple built in configuration tool allows the system to connect to almost any installed nurse call system.

Where complex systems are installed, a range of adaptors is available from your supplier.

When used in this way, the **Proxi-Mate Pro Nurse Call Set** supplements your existing nurse call pendant. The sensor activates the system just as if the resident had pressed the call bell themselves.



Proxi-Mate was designed and is manufactured in Australia and is fully compliant with all necessary standards.



SHA071114



Email your Purchase Orders to sales@statina.com.au or fax to 1300 308 259.

For more product information please phone Statina Healthcare on 1300 365 404.

www.statinahealthcare.com.au

Proxi-Mate



www.statinahealthcare.com.au

Proxi-Mate Bed and Chair Sensor

There are a number of different sensors available for the Proxi-Mate system.

The standard **Proxi-Mate Bed Mat Sensor**, for example, is commonly used when there is a need for a discreet and effective means of alerting staff when an individual has risen from their bed.

Likewise, the **Proxi-Mate Chair Mat Sensor** is designed to fit snugly on a chair squab.

Both of these sensors use state of the art, non-contact electronic proximity sensing. Unlike pressure sensors, these are completely reliable regardless of the weight of the individual. In addition, these sensors are so robust they come with a full 24 month replacement warranty. Many customers are still using their original sensors five years on and still going strong.



Proxi-Log USB Downloader

Both the **Proxi-Mate Pro Wireless** and **Nurse Call** systems record the date and time of each event as it occurs. The previous 100 events can be easily viewed directly on the devices' LCD screen and, while this is useful for a quick check of when things happened, it is difficult to form an overall picture when transcribing this data into a more readable format.

The **Proxi-Log USB Downloader** allows you to download the full data log (4000 events) to standard USB Flash drive. It consists of a USB adaptor combined with the **Proxi-Mate Report generator**.

The **Proxi-Mate Report generator** formats all the log data into a variety of pre-defined reports. This allows the user to scroll through each entry and to print documentation that can be used to support falls prevention or quality assurance activities. The data can also be used as part of a baseline assessment to determine when care staff intervention would be most effective.

Unlike traditional monitors, **Proxi-Mate sensors** are not 'single patient use' and so don't need to be replaced between patients/residents.

SHA071114