


# PRESSURE CARE

## Protect-A-Limb®



### Protect-A-Limb®

*The ultimate impact protection for frail skin*



Protect-a-Limb® protects the limbs from skin tears and skin damage. It is light in weight and yet bulky enough to protect. Aged Care Health Professionals highly recommend Protect-a-Limb® as a truly effective means of preventing recurring skin tears, with the resultant saving to healthcare costs. It is easy and fast to apply when compared to other methods of protection.

Protect-a-Limb is an Australian made healthcare product designed and tested for Australian conditions. It has been engineered to ensure the person requiring the garment remains cool in summer and warm in winter - a necessity in our harsh Australian climate.

The "cool in summer and warm in winter" element of Protect-a-Limb® is achieved by the unique Coolmax® System. Coolmax®, a registered trademark of INVISTA, is a channelled synthetic fibre on the inner of Protect-a-Limb®, surrounded by cotton on the outer surface which is slightly retentive.

#### Features:

- Light and comfortable
- Cool in summer, warm in winter
- Quick and easy to apply
- Dries faster than other fibres
- Is cost effective - saving time in application
- Gives protection against friction to vulnerable limbs
- Aids in patient comfort
- Coolmax® fabric wicks away moisture and prevents maceration
- Easy to clean - can be washed as easily as a pair of socks
- Finished to prevent fraying

Product Code	Product Name	Size	Length	Circumference	Use
PCMA-PAL-SM	Protect-A-Limb (Pair)	Small	40cm	18cm - 25cm	Arm & Small Leg
PCMA-PAL-SMBX12	Protect-A-Limb - Box of 12 Pairs	Small	40cm	18cm - 25cm	Arm & Small Leg
PCMA-PAL-MD	Protect-A-Limb (Pair)	Medium	50cm	25cm - 40cm	Arm & Average Lower Leg
PCMA-PAL-MDBX10	Protect-A-Limb - Box of 10 Pairs	Medium	50cm	25cm - 40cm	Arm & Average Lower Leg
PCMA-PAL-LG	Protect-A-Limb (Pair)	Large	60cm	32 - 44cm	Full leg
PCMA-PAL-LGBX8	Protect-A-Limb - Box of 8 Pairs	Large	60cm	32 - 44cm	Full leg

SHA010914-3